

³²P LABELED RIBOPROBES FOR NORTHERNS

1. Add the following in order:

H ₂ O	1.75 μl
1 mgml ⁻¹ BSA	2.5 μl
5 x buffer	5 μl
triton X-100	0.25 μl
1M DTT	0.5 μl
RNAguard	2 μl
10mM ATP	1 μl
10mM CTP	1 μl
10mM GTP	1 μl
100μM UTP	3 μl
DNA	1 μl (1 μg)
α ³² P-UTP 10mCi ml ⁻¹	5 μl
T3/T7/SP6 polymerase	1 μl

2. Incubate at 30°C for 3 hours.

3. Dilute into 100 μl by adding:

2 μl RNAGuard
10 μl 10 x DNAase buffer (manufacturer's)
10 μl 100mM DTT
52 μl H ₂ O

Remove 2 x 1 μl and spot onto DE81 paper for quantification later.

Add 1 μl DNAase and incubate for 15' at 37°C.

4. Phenol/chloroform extract.

5. Ethanol precipitate by adding 50 μl 7.5 M NH₄OAC, 2 μl 50 mgml⁻¹ tRNA and 376 μl ethanol.

6. Precipitate and resuspend in 500 μ l 50%formamide.

Quantification.

Rinse 1 filter: 5 x 5' in 0.5 M Na_2HPO_4
2 x 1' in dH_2O
2 x 1' in 95% ethanol

Dry thoroughly.

Count both filters in aqueous scintillant. Washed filter represents incorporated counts.