

## AFFINITY PURIFICATION OF ANTIBIODIES USING HIS-TAGGED PROTEIN

### THINGS YOU MUST KNOW FIRST:

- solubility of protein in TBS or Guanidine/Urea
- HIS-tagged protein does not elute from column in TBS or 4M MgCl<sub>2</sub>
- rough idea of capacity of Ni-agarose for antigen (normally 5-10mg per 1 ml resin). Want to aim to load resin to 80% capacity.

NOTE: It is possible to start this protocol with crude bacterial extract and essentially purify the HIS-tagged protein on the same column as you affinity purify the antibody. However, you get a better of idea of binding efficiency if you first purify the HIS-tagged protein, quantify the amount eluted and reload a new column as below.

1. Mix 8 mg purified HIS-tagged protein (that has been eluted from a Ni column in **Buffer E**) with **Buffer B**, to make a final volume of 8 ml. Adjust pH to 8.0 if necessary.
2. Add 2 ml 50% Ni-agarose slurry and shake at room temperature for 1h.
3. Decant protein/bead slurry into a column. Allow column to drain by gravity flow until just above top of beads.
4. Wash column with 10 ml (i.e. 10 column volumes) **equilibration buffer**.
5. Load 2 ml crude serum onto column. Collect flow through & reload 3 x. (Check activity of final FT on western - if it still reacts against the antigen, you have overloaded the column with serum).

6. Wash column with 5 ml (5 column volumes) **equilibration buffer** then 5 ml of: 2 M NaCl, 50 mM Tris pH 7.5.
7. Add 1 ml 4M MgCl<sub>2</sub>, allow to drain into column so that it has replaced previous solution then stopper column and leave for 15' at room temp.
8. Add more 4M MgCl<sub>2</sub> and collect 2 x 1 ml fractions as they drain from the column. Remove 10µl of each fraction for gel analysis.
9. Dialyse the combined fractions against dH<sub>2</sub>O for 1 h, then against PBS overnight.
10. Remove antibody solution from dialysis tube. Concentrate if necessary and store in PBS plus 0.2% azide.

NOTE: If you wish to reuse the column to purify more antibody against the same protein, wash the column with excess equilibration buffer. This should return the bead/protein complex to a state where it can be reused.

**Buffer B:**            8 M Urea  
                               0.1 M NaH<sub>2</sub>PO<sub>4</sub>  
                               0.01 M Tris-HCl, pH 8.0  
                               pH directly before use - 160µl 5M NaOH per 10ml

**Buffer E:**            As B, but pH 4.5  
                               pH directly before use - 11µl HCl per 10 ml

**Equilibration buffer:**

|                        |           |
|------------------------|-----------|
|                        | 100ml     |
| 150 mM NaCl            | 3 ml 5M   |
| 50 mM Tris-HCl, pH 7.4 | 2.5 ml 2M |