

ISOLATION OF MAIZE CHLOROPLASTS FOR PROTEIN IMPORT STUDIES.

Mark Settles (modified from Ken Cline)

1. Bring trays of seedlings from growth chamber to bench. Let sit in ambient light for ~1 hour. This reduces starch accumulation in the chloroplasts. The seedlings should be grown 10-14 days in soil without fertilizer in standard maize growth conditions (78°F day, 65°F night, 18 hour days).
2. If using continuous gradients, make a Percoll gradient by mixing 14 ml **2 x GR** with 14 ml Percoll in a SW28 Ultraclear centrifuge tube. Establish gradient by spinning at 17,000 rpm in SW28 (52,000 x g) for 40'. Keep the gradient cold (4°C in centrifuge, ice on bench).
3. The seedlings should have about three leaves showing. Cut the seedlings to get as many of the younger leaves as possible. Chop 25 g of leaf tissue into 200-250 ml **GR** (ice cold). You need about 40 seedlings for 25 g of tissue.
4. Grind the tissue in a blender with 4 very short pulses on a low speed. Turn the blender on then off almost immediately.
5. Filter through 1 layer of pre-wet miracloth or cheesecloth. Squeeze gently to increase the yield.
6. Spin down chloroplasts at 3,000 g for 8'. A swing out rotor is probably best for this step.
7. Resuspend the pellet in 5 ml of **GR** (more than this will overload a 28 ml gradient) using a wide bore transfer pipette or a camel hair brush. Layer onto the Percoll gradient.
8. Spin at 6,000 rpm (6,500 x g) in the SW28 for 15'.
9. Remove the top band (broken plastids) and discard. Take the lower band and dilute three fold with **import buffer (IB)**.

10. Spin down chloroplasts at 4,000 rpm (2,600 x g) for 8' in HB-4 rotor.
11. Resuspend the pellet in 15 ml **IB**.
12. Spin down the chloroplasts for 4' at the same speed as in 10.
13. Resuspend in an appropriate volume of **IB** (~1 ml of **IB** per 25 g tissue) to get ~1mg chlorophyll per ml.
14. Test chlorophyll concentration. Make a 1:10 dilution of the chloroplasts in **IB**. Mix the following in four tubes

Water	Chloroplasts	Acetone
200µl	0µl	800 µl
160µl	40µl	800 µl
140µl	60µl	800 µl
120µl	80µl	800 µl

Take the OD₆₅₂ of each tube. If the OD is above 1 dilute the chloroplasts more to get an accurate measure. Calculate chlorophyll concentration with the following formula:

$$\text{mg/ml} = (\text{OD}_{652})(0.02899)(1000/\text{ml of chloroplasts})(\text{dilution factor})$$

GR Buffer (1L)

50 mM Hepes-KOH pH7.5	50 ml 1 M Hepes
0.33 M Sorbitol	100 ml 3.3 M Sorbitol
1 mM MgCl ₂	1 ml 1 M MgCl ₂
1 mM MnCl ₂	1 ml 1 M MnCl ₂
2 mM EDTA	2 ml 0.5 M EDTA
5 mM Na-ascorbate	1 g Na-ascorbate
1% BSA	10 g BSA

Add the Na-ascorbate and BSA fresh the day you use the GR

Import Buffer (1L)

50 mM Hepes pH 8.0	50 ml 1 M Hepes
0.33 M Sorbitol	100 ml 3.3 M Sorbitol

If you are using step gradients, use an 85% step and a 35% step. For a 30 ml tube pour 6 ml 85% and 15 ml 35%. Increase these volumes if you are loading more than 3 ml on the gradient. The steps should be in GR.