

END LABELING DNA

1. Digest DNA in appropriate restriction buffer.
2. Add 1 μl shrimp alkaline phosphatase (SAP). Incubate **either** for 30' at 37°C if cut leaves 5' overhang **or** for 15' at 37°C and then 15' at 56°C if cut leaves 3' overhang.
3. Add 1 μl SAP and repeat incubation.
4. Add 2 μl 0.5 M EDTA. Incubate 10' at 68°C.
5. Extract twice with phenol/chloroform and then ethanol precipitate.
6. To 0.5 - 1.0 μg DNA add:
 - 5 μl **10 x kinase buffer**
 - 0.5 μl 10 mM spermidine
 - 25-50 μCi $^{32}\text{P}\gamma\text{ATP}$
 - up to 50 μl with dH_2O
7. Add 1 μl polynucleotide kinase and incubate for 5' at 37°C.
8. Add 2.5 μl 20 μM ATP and incubate for a further 30' at 37°C.
9. Ethanol precipitate to remove free nucleotides.

10 x kinase buffer (1 ml) (or use manufacturer's)

500 mM Tris pH 7.5	250 μl 2 M
100 mM MgCl_2	100 μl 1M
100 mM β -mercaptoethanol	6.9 μl