

FORMALDEHYDE RNA GELS

(Nelson protocol + modifications)

Always keep the RNA on ice.

1. Pretreat gel box, comb and tray with 0.1 M NaOH for 30' before use. **Rinse well with deionized water.**
2. Pour 1.5 % agarose gel in 1x Gel Buffer:
 - a) add 17 ml formaldehyde solution (37%) to 33 ml dH₂O in a flask in the fume hood. Mix.
 - b) rinse comb and gel tray well with deionized water and set up in the fume hood.
 - c) add 1.5 g agarose to 10 ml **10x Gel Buffer** and 40 ml dH₂O in flask. Microwave for 2'. Make sure agarose is completely dissolved and add to formaldehyde solution in fume hood. Pour gel immediately.
3. While gel is setting, rinse gel box well. Dilute **10 x Running Buffer** to 1 x and fill to the platform. The less buffer over the gel, the less diffusion of formaldehyde.
4. Prepare RNA samples by adding an appropriate amount of RNA and DEPC treated H₂O to a final volume of 7 μ l. Add 23 μ l **RNA Sample Buffer** to each RNA sample and mix well by pipetting up and down a few times. More or less RNA can be used but the 23:7 ratio must be maintained.
5. Heat samples to 65°C for 5' and place on ice. Load samples into the gel. The sample is heavy enough to sink in the well, but add tracking dye + loading buffer to the last lane to keep track of how far the gel has run.

6. Run the gel overnight. Place a glass plate on top of the gel to keep the formaldehyde from diffusing out and run at 25 volts. If running at high voltage (80-100v) keep the mA @ 75.
7. Photograph gel under UV and quantitate RNA loading using Biorad Multi-fluor S.

Solutions:

10 x Gel Buffer (100 ml)

200 mM MOPS pH 8.0	20 ml 1 M
50 mM NaOAC	5 ml 1 M
10 mM EDTA	2 ml 0.5 M

Final pH with added formaldehyde is 7. Autoclaving is OK but not necessary (the MOPS turns yellow when autoclaved, but it's fine). Wrap in foil.

10 x Running Buffer (2L)

200 mM MOPS pH 7.0	400 ml 1 M
50 mM NaOAC	100 ml 1 M
10 mM EDTA	40 ml 0.5 M

RNA Sample Buffer

300 μ l 10x Gel Buffer
1500 μ l formamide
14 μ l EtBr [10 mg ml^{-1}]
486 μ l formaldehyde

Mix in 15 ml tube and aliquot into eppendorf tubes. Store at -80°C and keep away from light. Use a 23:7 buffer to RNA ratio