

## DNA ISOLATION USING CTAB

1. Freeze plant tissue in liquid nitrogen and grind to fine powder.
2. Place a small spatula-full amount of ground tissue in 2 ml conical bottom eppendorf. Add 1 ml of **extraction buffer** (very viscous) and mix with yellow tip or cocktail stick.
3. Incubate at 65°C for 15 minutes. Mix by inversion after 5 minutes.
4. Add 714 µl of chloroform:isoamyl alcohol (24:1). Centrifuge at 8000 rpm for 5 minutes at room temperature
5. Transfer aqueous (top) layer to a 1.5 ml eppendorf and add 70 µl of 3M ammonium acetate. Mix by inversion.
6. Add 726 µl of isopropanol and mix by inversion.
7. Centrifuge at 7500 rpm for 5 minutes at room temperature, discard supernate and wash pellet in 600 µl of 70% ethanol.
8. Centrifuge at maximum speed for 15 minutes, discard ethanol and vacuum (or air) dry pellet.
9. Resuspend pellet in appropriate volume of 1 x TE buffer.

### EXTRACTION BUFFER (100 ml)

10% PEG 8000	10 ml
1 M Tris-HCl, pH 7.0	10 ml
5 M NaCl	28 ml
10% CTAB (heat to dissolve)	20 ml
500 mM EDTA	4 ml
10% SLS (sodium lauryl sarcosine)	3 ml
dH <sub>2</sub> O	25 ml