

IMMOBILIZATION ON CYANOGEN BROMIDE ACTIVATED SEPHAROSE

- Must have protein free of buffers that contain amino groups (that includes Tris) as these groups will attach to the resin and interfere with the coupling reaction.
 - 1g dry gel gives about 3.5ml swollen gel which can bind 30-200nmol ligand
 - want to aim for 80% binding capacity i.e. 160nmol. Sometimes extensive coupling will lead to loss of antibody binding sites. In which case, stop coupling reaction at 70% efficiency.
1. Gently suspend 1g CNBr-activated Sepharose in 1mM HCl (86 μ l concentrated in 1L) gently on a sintered glass funnel. Stir slurry with glass rod for about 20 minutes until matrix is swollen, then apply vacuum to remove liquid.
 2. Wash matrix with 3 x 100ml 1mM HCl. Wash to dryness each time.
 3. Wash matrix with 50ml **coupling buffer**.
 4. Dilute 30-200 nmol ligand in 15ml **coupling buffer**. Remove 100 μ l aliquot for later analysis.
 5. Quickly transfer activated Sepharose to a flask containing the ligand solution.
 6. Shake gently overnight at 4°C.
 7. Filter matrix on sintered glass funnel and collect flow through. Remove 10 μ l for gel analysis (compare to 10 μ l of starting solution to estimate binding efficiency). Keep remainder until you know whether the coupling reaction has worked.
 8. Wash matrix with 100ml **coupling buffer** on sintered glass filter.

9. Transfer matrix back to flask and incubate with 50 ml **blocking buffer** at room temperature for two hours or 4°C overnight.
10. On sintered glass filter, wash matrix with 50 ml **coupling buffer** and then 50 ml **acetate buffer**.
11. Wash with 4 cycles of 50ml **coupling buffer** followed by 50ml **acetate buffer**.
12. Finish the washes with coupling buffer containing 0.02% azide. The matrix is now ready for storage.

Coupling Buffer	100ml
0.1M NaHCO ₃	10ml 1M
0.5M NaCl	10ml 5M
pH8.3	

Blocking Buffer	100ml
1M ethanolamine	
pH8.0	

Acetate Buffer	100ml
0.1M sodium acetate	3.3ml 3M
0.5M NaCl	10ml 5M
pH4.0	