

IMMUNOLocalIZATION (FLUORESCCEIN)

1. Dewax slides 2 x 10' in histoclear.
2. Rehydrate by passing through the following ethanol series (2' each):
 - 99% ethanol
 - 99% ethanol
 - 95% ethanol
 - 85% ethanol
 - 70% ethanol
 - 50% ethanol
 - 30% ethanol
 - dH₂O
 - dH₂O
3. Incubate in **PBS/BSA** for 5' at room temperature.
4. Drain slides, place in a large petri dish containing wet tissues around the side. Put 50-100 μ l 0.1% Goat IgG in **PBS/BSA** on each slide. Cover petri dish and incubate for 15' at room temperature.
5. Drain slides and rinse in **PBS/BSA**.
6. Drain slides, place back in petri dish and put 50-100 μ l 1^o antibody on each (1/10 - 1/100 dilution in **PBS/BSA**). Incubate for 15' at room temperature.
7. Wash 2 x 15' in **PBS/BSA**.
8. Drain slides, place back in petri dish and put 50-100 μ l FITC conjugated 2^o antibody on each (1/10 - 1/100 dilution in **PBS/BSA**). Incubate for 15' at room temperature.
9. Wash 2 x 15' in **PBS/BSA**.
10. Stain with 100 μ M Evans Blue in **PBS/BSA** for 10' at room temperature. (optional).

11. Rinse in **PBS/BSA**.
12. Mount in 2:1 glycerol:PBS containing 2% N-propyl galate. Seal coverslips with nail varnish. View under fluorescent microscope with FITC filter.

PBS/BSA (1L)

| | |
|------------------------------|-----------|
| 0.15M NaCl | 30 ml 5 M |
| 10mM phosphate buffer pH 7.0 | 10 ml 1 M |
| 1mgml ⁻¹ BSA | 2 g |