

IMMUNOREACTION WITH WESTERN BLOTS

1. Remove filter from blotting apparatus and incubate in **blocking buffer** for 3 x 10' at room temperature with shaking.
2. Discard last wash and replace with antibody diluted in **blocking buffer**. Normally 10 μ l serum is added to 10 ml **blocking buffer**. Incubate shaking for 1 - 24 h. If incubating overnight, shake in cold room.
3. Pour off antibody into a falcon tube and store at -20°C. This solution can be re-used up to 20 times.
4. Wash filter in 3 changes of **blocking buffer** - 10' each.
5. Discard last wash and replace with 2° antibody diluted in **blocking buffer**. Normally 5 μ l 2° antibody is added to 10 ml buffer. Incubate at room temperature for 1 - 3 h.
6. Pour off antibody and wash filter in 3 changes of **blocking buffer** - 10' each.
7. Discard last wash and replace with **stain solution**.
8. When bands are developed, rinse filter in dH₂O to stop the reaction.
9. Filter can stored dry in saranwrap and used again at a later date.

Notes:

If using a secondary detection system that uses biotin you must use blocking buffer B. All other detection systems are OK with buffer A which is cheaper.

The stain solution used depends on whether horseradish peroxidase or alkaline phosphatase conjugated antibodies are used.

Blocking buffer A (2L)

5% non fat milk

100 g

| | |
|-------------------------------|-----------|
| 125 mM NaCl | 50 ml 5M |
| 10 mM phosphate buffer pH 7.0 | 20 ml 1 M |

Blocking buffer B (2L)

| | |
|-------------------------------|-----------|
| 1mgml ⁻¹ BSA | 2 g |
| 125 mM NaCl | 50 ml 5M |
| 10 mM phosphate buffer pH 7.0 | 20 ml 1 M |

Stain solution (HRP) (12 ml)

| | |
|--|--|
| 500 µgml ⁻¹ 4-chloro-1-naphthol | 2 ml of 3 mgml ⁻¹ in methanol |
| 125 mM NaCl | 300 µl 5 M |
| 10 mM phosphate buffer pH 7.0 | 120 µl 1 M |

Add 5 µl 30% H₂O₂ at the last minute.