

## INVERSE PCR (Susie)

1. Digest approximately 2  $\mu\text{g}$  of genomic DNA with the appropriate enzyme in a 50  $\mu\text{l}$  reaction. i.e.:

DNA	2 $\mu\text{g}$
10x buffer	5 $\mu\text{l}$
spermidine (100mM)	1.5 $\mu\text{l}$
restriction enzyme (10 U/ $\mu\text{l}$ )	2 $\mu\text{l}$
dH <sub>2</sub> O	to 50 $\mu\text{l}$

2. Clean the digest using a Wizard DNA Clean-up System kit (or other similar method - final volume should be about 50  $\mu\text{l}$ )

3. Self ligate 12.5  $\mu\text{l}$  of the cleaned DNA in a 100  $\mu\text{l}$  reaction. i.e.:

DNA	12.5 $\mu\text{l}$
10 x T4 DNA ligase buffer	10 $\mu\text{l}$
T4 DNA ligase	1 unit
dH <sub>2</sub> O	to 100 $\mu\text{l}$

4. Incubate at 12°C overnight.

5. The PCR is based on the Boehringer Mannheim Expand Long Template PCR System protocol. Set up two reactions, one with about 10  $\mu\text{l}$  of ligated DNA and the other with about 20  $\mu\text{l}$  of ligated DNA. A hot-start method is used, so use two 25  $\mu\text{l}$  mixes for each 50  $\mu\text{l}$  reaction, as follows:

### Mix 1

ligated DNA (see step 3)	10/20 $\mu\text{l}$
dNTPs (25 mM)	1.25 $\mu\text{l}$
each primer (10 $\mu\text{M}$ )	1.5 $\mu\text{l}$
dH <sub>2</sub> O	to 25 $\mu\text{l}$

### Mix 2

	10x buffer 3 (Expand Long Template System kit)	5
$\mu\text{l}$	Expand enzyme mix (DNA polymerases)	
0.75 $\mu\text{l}$		
	dH <sub>2</sub> O	19.25
$\mu\text{l}$		

6. Overlay Mix 1 with about 30  $\mu\text{l}$  paraffin oil, then heat the tubes to 94°C, in the first step of the PCR. Once the tubes are at 94°C, add Mix 2, beneath the oil. Then run the PCR:

35 rounds of    94°C 10 secs  
                      55°C 30 secs  
                      68°C 2 minutes  
 then                68°C 10 minutes

7. A nested reaction can then be carried out, either using 1  $\mu\text{l}$  of the first reaction as template, or by running the reaction on a gel, gel-extracting the band of interest and then using this as template. I found the second method gave a clearer, stronger band when the nested PCR products were subsequently run on a gel.