

PLATE METHOD FOR λ DNA PREP

1. Plate 4 plates to confluence. Incubate 37°C, overnight.
2. Add 5 ml of SM per plate. Swirl plates gently for 2 hrs at room temperature.
3. Pool the liquid from the plates in a 50 ml tube. Add 500 μ l of chloroform (crude cleaning step), vortex, and leave for 30 minutes (vortexing every now and then).
4. Transfer contents to a fresh tube, leaving chloroform behind. Add 6 μ l of 10 mgml⁻¹ RNAase and 4 μ l of 10 mgml⁻¹ DNAase. Leave at room temperature for 30 minutes.
5. Precipitate the phage by adding 1.2 g of NaCl. Dissolve salt and leave at 4°C for at least 1 hr.
6. Transfer sample to 30 ml tube and spin in SS34 rotor at 7000 rpm, 4°C, for 20 minutes.
7. Transfer supernatant to a fresh 30 ml tube and add 2 g PEG6000. Dissolve the PEG (takes 10-30 minutes to fully dissolve), and leave sample at 4°C overnight.
8. Spin at 7000 rpm for 20 minutes, as before. Remove as much supernate as possible and discard. Resuspend pellet in 500 μ l of SM and transfer to eppendorf.
9. Extract the DNA, using first an equal volume of chloroform, then phenol, then phenol-chloroform, and lastly chloroform again.

10. Precipitate DNA with 2 volumes of 100% EtOH, spin and wash with 70% EtOH. Dry pellet and resuspend in appropriate volume (about 100 μ l) of 1x TE buffer.