

MINI PHAGE PREPS.

1. Inoculate 5 ml LB with 0.2 ml plating cells and 300 μ l phage stock. Incubate at 37°C overnight.
2. Pour 1.3 ml phage culture into an eppendorf tube, microfuge 1' and transfer supernate to a fresh tube which contains 300 μ l of 20% PEG 8000, 2.5 M NaCl. Mix and leave at 4°C for more than an hour.
3. Microfuge for 5' at 4°C. Resuspend pellet in 100 μ l 300 mM NaCl, 100 mM Tris pH8, 1 mM EDTA.
4. Add 200 μ l phenol:chloroform 1:1, mix and sit for 5'. Spin in microfuge and transfer aqueous layer to a fresh tube.
5. Add 10 μ l 4.4 M ammonium acetate pH 5.2 and 400 μ l ethanol.
6. Mix and microfuge for 5'.
7. Wash pellet in 70% ethanol, dry and resuspend in 20 μ l TE.
8. Use 10 - 20 μ l for a restriction digest. (Adding RNAase to the digest).