

PLASMID DNA MINIPREP

1. Use 1 ml of an overnight liquid culture or a colony scrape off a plate.
2. If using a liquid culture, pellet cells in a microfuge.
3. Add 100 μ l **Solution I**. Stand at room temperature for 5'.
4. Add 200 μ l **Solution II**, shake gently and place on ice for 5 - 10'.
5. Add 150 μ l **Solution III**, vortex and place on ice for 10'.
6. Spin down and remove supernate to fresh tube.
7. Add 750 μ l ethanol to the supernate, incubate at room temperature for 5' and then pellet.
8. Resuspend in TE, add 1 μ l 10 mgml⁻¹ RNAase A, incubate for 10' at 37°C.
9. Phenol/chloroform extract, ethanol precipitate and resuspend at concentration required for use.

Solution I (100ml)

10 mM Tris pH 8	500 μ l 2 M
1 mM EDTA	200 μ l 0.5 M

Solution II (50 ml) (Make fresh each time)

0.2 M NaOH	10 ml 1 M
1% SDS	5 ml 10%

Solution III (100 ml)

3 M K Ac pH4.8	29.4 g
----------------	--------

(dissolve in as little water as possible. Needs about 12 ml acetic acid to pH)