

PLATING & TITERING λ PHAGE

Preparation of plating cells.

1. Streak a colony of host cells on LB and incubate at 37°C overnight.
2. Inoculate 50 ml LB + 10 mM MgSO₄ + 2% maltose, incubate 4-6 h at 37°C or overnight at 30°C. Do not let OD₆₀₀ exceed 1.
3. Pellet bacteria at 2000 rpm for 10'. Resuspend cells in half the original volume of sterile 10 mM MgSO₄.
4. Dilute cells to an OD of 0.5 with 10 mM MgSO₄. Store at 4°C and use within 48 h.

This is the correct way to do it, however, you get just as good results using an overnight culture that has been grown as in 2 and then stored at 4°C.

Plating.

1. Make dilutions of phage stock in SM.
2. Mix 200 μ l plating cells with up to 100 μ l phage solution. Incubate without shaking at 37°C for 15'.
3. Add 5 ml top agarose (LB + 0.7% agarose), mix well and pour over an LB plate. Incubate 37°C overnight.
4. Calculate titer on the basis of pfu (plaque forming units) per μ l phage stock.

SM (100 ml)

2 ml 5 M NaCl

2.5 ml 2 M Tris pH 7.5

1 ml 1 M MgCl₂

500 μ l 2% gelatin