

RANDOM PRIMING DNA

(Feinberg and Vogelstein (1983) Anal. Biochem. 137, 266-267.)

Can use purified fragment in which case start at 5.

1. Digest DNA in appropriate restriction buffer. Electrophorese on a 1-2% low melting agarose gel in cold room.
2. Excise fragment of interest and weigh agarose slice in a preweighed eppendorf tube (trim as much agarose as possible).
3. Add H₂O at a ratio of 3 ml H₂O to 1 g agarose.
4. Boil tube for 7' and then measure volume of liquid. Calculate concentration of fragment per ml of solution based on amount of DNA fragment in gel (assume 100% recovery). Store at -20°C in aliquots of 50 µl.
5. Use approximately 20 ng DNA per labeling reaction. Boil DNA for 3' and place on ice if purified fragment or at 37°C if in agarose.
6. Set up labeling reaction as follows:

15.5 µl DNA and H₂O
5 µl **OLB**
1 µl 10 mgml⁻¹ BSA
1 µl Klenow polymerase
2.5 µCi ³²PαdCTP
7. Incubate at 37°C for 1- 3 hours or overnight at room temperature.
8. Add 300 µl TE and boil for 10' prior to use.

To make **OLB**:

Solution A (Store at -20°C)

625 μ l 2M Tris-HCl pH 8
125 μ l 1M MgCl₂
18 μ l β -mercaptotethanol
5 μ l 100 mM dATP in TE
5 μ l 100 mM dTTP in TE
5 μ l 100 mM dGTP in TE
250 μ l dH₂O

Solution B (Store at 4°C)
2 M HEPES pH 6.6 with NaOH

Solution C

Hexadeoxyribonucleotides (Pharmacia #2166) evenly suspended in TE at 90 OD unitsml⁻¹.

OLB = **A** 100 : **B** 250 : **C** 150

Store at -20°C in 50 μ l aliquots.