

# SOUTHERN BLOTTING

1. Run DNA on gel.
2. Photograph gel.
3. Denature DNA by shaking gel in 1.5 M NaCl, 0.5 M NaOH for 45'-1 hr.
4. Neutralize gel by shaking in 3 M NaCl, 0.5 M Tris pH7.5 for 45' - 1 hr.
5. Cut a piece of Nytran and two pieces of whatman filter paper the same size as the gel and wet in 2 x SSC.
6. Set up blotting tray with **20 x SSC** in tray as follows:



7. Leave overnight at room temperature.
8. Remove filter and rinse for 5' in 2 x SSC.
9. X-link filter in Stratalinker.
10. Either use immediately or store at room temperature between whatman filter paper.

## 20 x SSC (2L)

3 M NaCl	351 g
0.3M sodium citrate	176 g
pH to 7.6 with HCl	

