

HYBRIDIZATION OF SOUTHERN BLOTS.

1. Prehybridize filter for 1 - 2 hr at 65°C in **prehybridization buffer**.
2. Hybridize overnight at 65°C in **hybridization buffer**.
3. Wash filters at 65°C for 2 x 30' in **1 x wash buffer** and if necessary for a further wash in **0.1 x wash buffer**.
4. Air dry filter and expose to film.

Prehybridization buffer (100ml)

3 x SSC	15 ml 20 x
1% SDS	10 ml 10%
0.1% sodium pyrophosphate	1 ml 10%
5 x Denhardt's	5 ml 100 x
200 µgml ⁻¹ sheared ss DNA	2 ml 10mgml ⁻¹
	67 ml H ₂ O

Hybridization buffer (100ml) (Don't use dextran for non-genomic blots)

3 x SSC	15 ml 20 x
1% SDS	10 ml 10%
0.1% sodium pyrophosphate	1 ml 10%
5 x Denhardt's	5 ml 100 x
200 µgml ⁻¹ sheared ss DNA	2 ml 10mgml ⁻¹
10% dextran sulphate	10 g
	up to 100 ml with H ₂ O

1 x wash buffer (2L)

1 x SSC	100 ml 20 x
1% SDS	20 g
0.1% sodium pyrophosphate	20 ml 10%

0.1 x wash buffer (2L)

0.1 x SSC	10 ml 20 x
-----------	------------

1% SDS	20 g
0.1% sodium pyrophosphate	20 ml 10%

100 × Denhardt's

10g Ficoll
10 g polyvinylpyrrolidone
10 g BSA
up to 500 ml with H₂O

Sheared salmon sperm DNA

Dissolve herring sperm or salmon sperm DNA at 10mgml⁻¹ in dH₂O.
Autoclave for 5' and rapidly cool. Store frozen in aliquots.

RECYCLING MEMBRANES

1. Place filter in 0.2M NaOH, 0.1%SDS for 30' at 37°C.
2. Rinse in 2 × SSC and then prehybridize as normal.