

WESTERN BLOTTING

1. Equilibrate gel in **blotting buffer** for 20' at room temperature.

2. Make sandwich as follows:

grey panel

sponge

whatman paper

gel

0.2 μ m nitrocellulose (or nylon)

whatman paper

sponge

clear panel

Make sure no bubbles are present in any of the layers.

3. Place blotting sandwich in the tank with the grey panel next to the cathode (black).

4. Fill tank with **blotting buffer**, add ice pack and stirbar.

5. Run at 30V; 40 mA overnight or 100V; 250 mA for 1h. Stir buffer throughout run.

6. Disassemble sandwich and either store blot between filter paper at room temperature or use immediately.

Blotting buffer (1L)

25 mM Tris pH 8.3 3.03 g Tris

192 mM glycine 14.4 g glycine

20% methanol 200 ml methanol